Ref.: 400060004

# Almohadilla eléctrica cervical

## Manual de instrucciones

Instructions manual  
Notice d'utilisation  
Manual de instruções

## Almohadilla eléctrica cervical

Cervical electric heating pad  
Coussin chauffant cervical électrique  
Almofada elétrica cervical

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## How to use

### Before first use

1. Remove all packaging.
2. Connect the controller to the heating pad by inserting the plug in coupler.
3. Plug the mains plug into the wall socket and switch on.

### Using your Pad

Before turning on the heating pad, make sure the switch is in the 0 position and plug it into the power outlet.

To turn on the heating pad, adjust the sliding switch to one of the temperature settings:

- 0: Switched off.
- 1: Lowest heat level.
- 2: Middle heat level.
- 3: Highest heat level.

[Note: the fastest way to heat up the pad is to switch it straight to number 3]

After a few minutes, the pad will have reached the desired temperature.

Place the heating pad using the strap and Velcro for the neck.

Connect the plug to the socket at the bottom of the heated vest and the other end to a power supply (AC220-240V). Ensure there is an electrical outlet nearby and that the power cord is not in a position where the user might trip over it.

After preheating, turn the control knob to the heat setting you need.

Place it on the part of the body you want and select the ideal temperature. Never use it with animals.

For enhanced heat effect, keep the pad covered to prevent heat from dispersing during intervals when the device is off.

To turn off the device, turn the knob to the "Off" position.

The heating pad can be used for up to 90 minutes before it automatically shuts off. If you wish to extend this time, it is advisable to let the pad cool for 5 minutes before using it again.

Once you have finished using the heating pad, slide the switch to position 0. Allow it to cool completely before storing.

### Periodic Checks

When the appliance is not in use, it is advisable to store it in its original packaging in a dry place and avoid placing heavy objects on it. Prevent folds from forming. Frequently examine the device and the flexible cable to ensure they do not show signs

of wear or damage. If any signs of wear or malfunction occur, take the appliance to any authorized service center.

## Cleaning

Small stains can be cleaned with a cloth or a damp sponge and, if necessary, with a small amount of liquid detergent for delicate fabrics.

### Warning

The electric pad should not be dried in a dryer or ironed. If the electric pad is very dirty, it can be hand washed or machine washed by removing the detachable controller and setting the machine to the 30°C program for delicate fabrics.

Use a neutral or delicate fabric detergent according to the manufacturer's instructions. Do not wash the electric pad too frequently to avoid excessive wear. The electric pad should not be washed in the washing machine more than five times during its lifespan.

When drying the electric pad, do not hang it with clothes pegs or other similar devices. Do not expose the pad to heat sources such as the sun or a radiator.

Reconnect the power cable to the electric pad only when both the plug and the pad are completely dry. Never switch on the pad to dry it.

## Notes

### Important

This pad is equipped with a safety device that automatically regulates the temperature by connecting and disconnecting the power at regular intervals. The temperature of the pad is the maximum allowed by current safety standards to ensure user safety. The device provides gentle heat that is not harmful.

## Malfunctions & Repair

1. In case of a malfunction, take the device to an authorized service center.
2. Do not attempt to disassemble or repair the device unless you are a qualified professional.
3. If the power cord is damaged, it must be replaced immediately.

## Caution!

This device can be used by children aged 8 and older, and by persons with physical, sensory, or mental impairments, or lack of experience and knowledge, provided they have been supervised or instructed on how to use the device safely and understand the associated risks. Children should not play with the device. Cleaning and maintenance should not be carried out by children without supervision.

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### Safety precautions

- 1.** Do not use by persons insensitive to heat and other very vulnerable persons.
- 2.** Do not use by children under the age of three (from 0 to 3 years).
- 3.** This electric heating pad is not to be used by young children (from 3 to 8 years) unless the controls have been pre-set by a parent or guardian or unless the child has been adequately instructed on how to operate the controls safely.

**4.** This electric heating pad can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.

**5.** Children shall not play with the electric heating pad.

**6.** Cleaning and user maintenance shall be made by children without supervision.

**7.** After removing the wrapping, make sure that the appliance is complete without visible damage that

could have been caused during transport. If in doubt, do not use the device. Contact the authorized technical assistance services.

**8.** Before connecting the device, make sure that the data on the board corresponds to the data on the electrical network. The plate is located in the device.

**9.** Do not use adapters, multiple plugs and / or extension cords.

**10.** This device must be used exclusively as heating pad and always following the instructions for use. Any other use will be considered inappropriate and, therefore, dangerous. The manufacturer declines

all responsibility for possible damages arising from improper and incorrect use.

**11.** The use of any electrical device entails compliance with some fundamental rules. Especially the following:

**a.** Never touch the appliance with wet or damp hands.

**b.** Do not pull on the power cord or the appliance to unplug it from the socket current.

**c.** Do not leave the device exposed to atmospheric agents such as rain or sun.

**12.** Before carrying out any cleaning or maintenance operation, disconnect the appliance from the mains, unplugging it.

**13.** In the event of failure and / or malfunction of the device, turn off and do not manipulate. For possible repairs, contact the authorized technical assistance services.

**14.** To avoid dangerous overheating, it is recommended to completely unwind the power cord and unplug the appliance from the mains when not in use.

**15.** Never change the power cord of this appliance. If it is damaged, the appliance should be thrown away.

**16.** When you decide not to use the appliance anymore, it is recommended

to disable it by cutting the power cable, after unplugging it from the power supply. Also, it is recommended to leave unused those parts of the device that may be dangerous, especially for children who could use the device out of use for their games.

**17.** Never immerse the device in water.

**18.** Prolonged application of the appliance at maximum temperature can cause skin burns. Do not fall asleep with the pad running. Do not pull the device by the power cord and do not use the cable as a handle. This device can not be used in

hospitals. Do not use the pad with patients, children or people sensitive to heat.

**19.** Do not use on very young children or people insensitive to heat.

**20.** Do not clean or dry clean.

**21.** Not suitable for commercial purposes.

**22.** The heated pad will automatically switch off if there is a fault.

**23.** The heated pad is fitted with a 90 minute timer.

**24.** The heated pad contains no user serviceable parts. If the heated pad stops working please contact a qualified repair person.

**25.** Electric and magnetic fields from this electric heating pad

may interfere with the operation of pacemakers. However, they are well below the limits: power of electric fields: 5 000 V/m max., magnetic fields power: 80 A/m max., magnetic induction: 0.1 mT Max please therefore consult your doctor and the manufacturer of your pacemaker before using this electric heating pad.

**26.** The electronic components in the switch will become hot when the heating pad using.

**27.** You should therefore not cover switch or place it on the heating pad when the latter is used.

**28.** This electric hea-

ting pad must be regularly checked for the presence of any sign of wear or damage. Where appropriate, this heating pad has been used inappropriately or if it heats more, it must be checked by the manufacturer before being used again.

**29.** If the heating pad power supply cable is damaged, it must be replaced by the manufacturer, its service agent or similarly qualified persons in order to avoid hazards.

**30.** If you want to store or fold the heating pad:

- Allow it to cool down before folding
- Do not crease it by placing items on the top of it during storage.

**31.** Do not use on adjustable bed or check that the heating pad or cord do not become trapped or rucked.

**32.** When using the electric heating pad for several hours, it is recommended to set the switch to the lowest temperature level to avoid overheating.

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---

#### **CORRECT DISPOSAL OF THIS PRODUCT.**

This marking indicates that this product should not be disposed with other household wastes throughout the EU.

To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmental safe recycling.

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**EN.** This product has a 3-year warranty from the date of sale, declining all responsibility for defects or breakages caused by misuse. For this guarantee to be valid, it is essential to present this card as well as the purchase receipt or invoice.

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